



EM-CONTOURING

A simple and fast way to build muscle and burn fat!




30 = 20,000


Mins Treatment Workouts

16% ↑ Average increase in muscle mass 19% ↓ Average fat reduction

BODY CON-TOURING MUSCLE



NON-INVASIVE



NON-RADIATING



NON-THERMAL



BODY SCULPT



NON-INVASIVE
BUTTOCK PROCEDURE



NO DOWNTIME



NO PAIN



OPERATOR
INDEPENDENT



NO CONSUMABLES



E-mail:
info@kmss.se



Tel:
+46 70-727 77 58



Web:
www.kmss.se



Application

Muscle Contractions

Female who seeks improvement in muscle and fat----BUTT LIFT, abs.

Models, Idols, Celebrities ---- contouring a sexy figure




male who seeks improvement in muscle and fat---muscle building, Chocolate abs.







Theory



EM-Contouring with HIFEM(High-Intensity Focused Electromagnetic) technology induces supramaximal muscle contractions not achievable by voluntary action. When exposed to supramaximal contractions, the muscle tissue is forced to adapt to such extreme condition.



Voluntary contractions may increase the demand for energy from fat cells. Released epinephrine signals the fat cells to initiate lipolysis. The fat stored in form of triglycerides is decomposed into free fatty acids (FFAs) and glycerol which are used as energy sources.

Intensive supramaximal contractions enhance release of epinephrine which triggers a cascade effect leading to supramaximal lipolysis in fat cells. Because the stimulation is so rapid and intensive, the lipolytic response is excessive and FFAs start over-accumulating in adipocytes.

The FFA overflow causes cell dysfunction and induction of apoptosis – programmed cell death. Dead cells collapse and are naturally flushed away.



Clinic Data

A course of treatment is 4 times,
each time only takes 30 minutes.
Do it at least 2 times a week and
2 weeks in a row easy and fast.

Average effect of 22 candidates:

18.6%

The fat thickness decreases

15.4%

The abdominal muscle thickness increases

10.4%

The Diastasis rectus
abdominis(DRA) improves

7%

The number of muscle
bundles increases

3.8cm

The waist
circumference lost



Female | Age 42
BIM25.1KG/M2(before)
BIM 24.4kg/M2 (after 2 months)
Weight -4.4kg 1b (-2.9%)
Fat thickness -32.9%
Muscle thickness +19.4%
DRA -15.9% waist circumference -5.7cm



Before



After



Before




After



Technical Advantages

01 Scientific Magnetic Pulses



02 Optimized Applicator

- Latest generation technology
- High-tech nonwater cooling
- High quality insulated material
- One size suitable for multiple areas



03 Automatic Control System

- Auto-adjustment energy range when changing frequency
- Smart control of total Tesla per session



04 User-Friendly Screen

- High Resolution Interface
- Free Sliding operation
- 45° adjustable for best viewing



05 4 Working Modes for simple use and professional use

- Simple 1 - For people who seldom exercise with less muscle and more fatness
- Simple 2 - For people who do regular exercise with more muscle and less fatness
- Specialty 1- For professional practitioners to customize protocol A
- Specialty 2 - For professional practitioners to customize protocol B

Specification

Screen	15” Sliding Screen
Pulse Frequency	1-150Hz
Energy Range	0-100%
Pulse Width	285us
Magnetic Intensity	0.4-2.5T
Treatment Time	Up to 60mins
Net Weight	31kgs

Benefits



Abdomen



Buttocks



Legs



Biceps beachii



Triceps brachill



E-mail:
info@kmss.se



Tel:
+46 70-727 77 58



Web:
www.kmss.se