K.M.S.S

Kanaan Medical Supplies & Services



CooL 360°

Freeze Away Stubborn Fats

E-mail: info@kmss.se

Tel: +46 70-727 77 58

Web: www.kmss.se

Cool 360° technology

New Cool 360° surrouding cooling system gives even cooling 60% faster from all angles, improving efficiency and results.

The larger encompassing of target area, more fat cells can be targeted and killed per treatment.









Faster, 35 minutes saved

Increased comfort

60% more tissue treated


E-mail: info@kmss.se
Tel: +46 70-727 77 58
Web: www.kmss.se


Theory

Cryolipolysis is

Harvard University patent technology, frozen de-fatting lock temperature technology, getting -10 to +5 °C under the skin, fat cells aging ahead of time and death. The fat cells in the treated area are gently dissolved through the body's normal metabolism process and then the dead cells in 1-3 months via lymphatic system.

0 DAY **30** DAYS **90** DAYS





0 Day 15 Days 30 Days 60 Days 90 Days

Many of us have stubborn fat despite diet & exercise.

Cryolipolysis technology uses controlled cooling to target & kill only these fat cells.

In the weeks to follow, your body naturally processes the fat and eliminates these dead.

Cryolipolysis procedure results are long-term, as treated fat cells are gone for good.

E-mail: info@kmss.se
Tel: +46 70-727 77 58
Web: www.kmss.se

Application

Treatment



Abdomen



Double Chin



Upper arms



Thighs



Buttocks



Waist

Clinical






Female | Age 36
Weight: 58.2kgs
Circumference: 56cm(left) 55.5cm(right)
(before)

inner thigh 2 treatment areas and 3 sessions
Weight: 54.1kgs
Circumference: 53.6cm(left) 51.8cm(right)
(after)

Female | Age 42
Weight: 70.3kgs
Circumference: 97.6cm(before)

Waist :4 sessions
Weight: 60.6kgs
Circumference: 82.2cm (after)





E-mail: info@kmss.se
Tel: +46 70-727 77 58
Web: www.kmss.se

Technical Advantages

01 Proven technology


Fat Reduction VS Weight Loss

They may sound similar, but fat reduction and weight loss are two very different cocept.

Losing weight doesn't always lead to a more sculpted appearance—but fat reduction can reshape and contour your body.

Fat Cells(After Weight Loss)
Fat cells shrink in volume

Fat Cells (After Cryolipolysis)
Fat cells numbers reduction



- What happens when you lose weight?


When you lose weight, the size of your fat cells decreases but the number of fat cells stays the same.


- What happens after your cryolipolysis treatment?

After treatment, targeted fat cells are eliminated for good so you will have fewer fat cells in the treated areas. Untreated areas will have no change in fat cell distribution.

E-mail: info@kmss.se
Tel: +46 70-727 77 58
Web: www.kmss.se

02 Cover more tissue







NEW Innovation 360 touched cooling.
100% effective cooling area, surrouding cooling action in entire applicator

V S

Traditional Cryolipolysis by Two Chilling Plates Only. 40% effective cooling area, cooling by two pads inside applicator





03 Localized fat tissue treatment ensure maximum safety.





360°surround cooling non invasive cold therapy fat reduction procedure makes fat cells under a process “apoptosis”, the applicator localized reduction of fat deposits. The degree of exposure to cooling causes cell death of subcutaneous fat tissue, without apparent damage to the overlying skin.

E-mail: info@kmss.se
Tel: +46 70-727 77 58
Web: www.kmss.se


04 Customized Applicators

HANDLE A

(Large Size)

24.6*62*48cm






HANDLE B

(Medium Size)

20.6*62*48cm






HANDLE C

(Small Size)

17*62*48cm





HANDLE D


(For double chin)


8.4*26*25cm

E-mail: info@kmss.se
Tel: +46 70-727 77 58
Web: www.kmss.se

06 User friendly software

- Intelligent parameter recommending.
- Each handle，after choosing gender and treatment area, the parameter will be recomended automatically.






- 4 handle independent pressure control
- Each handle, the cooling and pressure control seperately, more comfortable and safer

Specifications

Power Supply Output	3000W
Handles Number	4
Vacuum Pressure	0-60kpa
treatment time	0-60mins adjustable
Cooling System	water+wind+semiconductor
Cooling Temperature	-10℃to10℃
Screen	15”true color touch LCD display
Language	multi languages
Voltage	110V/220V
Package size	67*63*123CM 76*51*44CM



Accurate Temperature Control System

New Technology 360 Degree Cooling

10 Size Contours Versatile use

Non Invasive No Pain No Rebound

Up to 40% fat reduction Satisfied Results

E-mail: info@kmss.se
Tel: +46 70-727 77 58
Web: www.kmss.se